Art. 2432 - 2433

Panacea "la petite" Electronic

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DRY STERILIZER- USER MANUAL

Remarks:

- ◆ These operator's instructions have to be carefully kept near the equipment, in order to be consulted by all users.
 - ◆ The use of this device is limited to qualified staff
 - ◆ Before the start, please check that your load withstands without any damages at the set temperature (+10% for safety).
 - ◆ The equipment is to be used dry, do not introduce water or other liquids.
 - ◆ In case of black-out the sterilizing cycle is not valid
 - ◆ In case of lack of voltage or lighting of the alarm pilot lamp the equipment stops the sterilization cycle.
 - ◆ We recommend not to lean any object on the equipment
 - ◆ The material should be placed in such a way that the hot air can flow freely around the material that is to be sterilised.
 - ◆ We recommend to clean and dry carefully the material to sterilize.
 - ◆ Temperature oscillations of few °C, during the cycle, have to be considered normal.
 - ◆ In case more consecutive cycles the apparatus could achieve elevated external temperatures.
- The equipment is only designed for the sterilisation of metal objects (scissors, stitch-removers, nail clippers etc), glass or other material capable of withstanding the temperature of the machine.
- ◆ It is inadvisable to place gauze, textiles and packaging inside for sterilisation.

Description of the equipment :

Referring to the figure of page 2 the sterilizer components are described here below:

- (A) **Master switch**
- (B) **"Ready machine" yellow pilot lamp**
This lamp is on when the sterilizer is ready and is waiting the start of a cycle
- (C) **Alarm red lamp** 
This pilot lamp is turned on in order to signal some possible anomalies (for ex. lack of voltage). When the master switch is on, the pilot lamp is turned on to signal the previous switching off. By pressing the START key the alarm is cancelled. If, during the cycle, this pilot lamp is turned on, the cycle is to be considered null (see also Alarms Table).
- (D) **START push button**
The function of this push button is to start a cycle, if the pilot lamp (B) is turned on.
The **START** push button is also working as RESET key, if it is pressed during a sterilization cycle the cycle is cancelled. Moreover, if an alarm is intervening, by pressing the **START** push button, the alarm is cancelled.
- (E) **Switch push button** 
The function of this push button is to alternate, on the display (H), the visualization of temperature and sterilization time.
- (F) **Temperature pilot lamp** 
When this pilot lamp is turned on it shows that the value read and/or set on the display (H) is referred to the temperature.
- (G) **time pilot lamp 1/2**
When this pilot lamp is turned on it shows that the value read and/or set on the display (H) is referred to the sterilization length of time, heating excluded.
- (H) **display**
Cycle parameters are visualized on the display : temperature inside the chamber, sterilization time, set temperature.
In order to reset the cycle data, this must be stopped.
 - The temperature can be set between 140°C and 200°C,
If the temperature inside the chamber is lower than 120°C you can read "t.Lo", (low temperature) on the display, if the temperature is higher than 220°C you can read "t.Hi" (high temperature).
 - Time can be set between 1 minute and 120 minutes, or as infinite time.
Time can be set between 1 minute and 120 minutes, it is moreover possible to set an infinite time that is anyway not recommended for high sterilization temperatures. In order to set an infinite time it is necessary to go beyond 120 minutes on the display. When the following symbols - - - will appear on the display wait 5 seconds for automatic storage.

ATTENTION: AT INFINITE TIME THE EQUIPMENT CAN BE SWITCHED OFF ONLY MANUALLY.

- (I) (L) **keys**  
 The function of these keys is to modify the temperature setting and the sterilization time, increasing or decreasing the value. By pressing the **+** key once, with the sterilizer waiting (B pilot lamp turned on) the set temperature value or the set sterilization time is visualized on the display.
 From this moment, by pressing the **+** or **-** keys it is possible to modify the value.
 After reaching the desired value, wait 5 seconds for automatic acquisition.
 If one of these keys is pressed during the cycle, the set values of time or temperature can be read on the display but it is not possible to modify them.
- (M) **Heating pilot lamp**
 This lamp is on when resistances are on. After the starting of a cycle, resistances are warm until the set temperature is reached. When the temperature is reached the **HEATING** pilot lamp is turned off and it is alternately turned on, as the equipment is keeping the sterilization temperature .

User's instructions

1. Place the material to be treated inside the sterilizer and close the door (do not open the door during the cycle).
2. Connect the plug into the 230V outlet. Turn the master switch on (A).
3. The "alarm" red pilot lamp is turned on to show the lack of voltage (given by the previous switching off of the machine).
4. Press the "START" push button once, in order to cancel the alarm.
 The yellow pilot lamp (B) is turned on at this moment, in order to show that the machine is ready for the cycle.
5. Set desired temperature and time of sterilization,
 Set the temperature : with the temperature pilot lamp on, press the **+** or **-** keys, or
 Keeping the push button pressed, the increase or the decrease of the set value becomes quick. After reaching the desired value, wait 5 seconds for the acquisition of the set value.
 Set the sterilization time:
 press the switch push button **6** turning on the timer pilot lamp $\frac{1}{2}$. Repeat the previous operations, in order to visualize and modify the time. Wait 5 seconds for the acquisition of the set value.
6. After the parameters setting, start the sterilization cycle by pressing the "START" push button. The yellow pilot lamp (B) is turned off and, if the temperature inside the chamber is lower than the set temperature, the "heating" pilot lamp (M) showing the resistances heating is turned on.
7. The sterilizer arrives at the set temperature and keeps it for the set time. The temperature measured inside the sterilization chamber is visualized on the display (with a + 10°C tolerance) or, by switching with the **6** push button, the remaining time to the end of the cycle.
8. At the end of the set time, the cycle automatically stops and the display shows the writing END, informing that the sterilizer has finished the cycle.

Important: the cycle is considered valid only if, at the end, the writing END is appearing.

9. If the equipment is not used, turn the master switch (A) off.

** Sterilization cycles

Temperature	Time in minutes
180°C	*45'
170°C	*75' (1h 15')
160°C	*120'

* minimum recommended time


ALARM TABLE

Display/Defect	Cause	Remedy
"T.Lo" / inside temperature too low. The minimum heating parameters have not been reached.	1) lack of voltage	1) repeat the cycle
	2) excessive load	2) reduce the load and repeat the cycle
	3) resistances breaking	3) replace the faulty resistance and repeat the cycle
"T.Hi"/ High temperature beyond 240°C	Thermostat breaking	1) open immediately the door 2) replace the thermostat
The equipment is completely stopped	1) Over-temperature, the safety thermostat is working	Open immediately the door 1) Reset the safety thermostat 2) Replace thermostat
	2) The fuses is broken	Replace the fuses
"CnF"	Faulty card	Replace the electronic card

(***) **SAFETY:** Should the safety thermostat intervene, all power goes off except to the green main switch. In order to reset the sterilizer it is necessary to unscrew the black cap situated on the back side of the machine and push the button its full course.

The overheating cut-off is operating because of an excess temperature in the chamber.

If during the following cycle the overheating cut-off should work again, call the technical assistance in order to eliminate the inconvenience causing its intervention.

(***) if existing in the version

Cleaning

To carry out a good cleaning use only neutral detergents and rinse with plenty of water.

The gaskets and all plastic parts can be cleaned with a cotton flock soaked in alcohol.

Do not use sponges or abrasive materials which can damage the surface.

To clean do not use detergents, disinfectants or other liquids containing chlorides as they are corrosive for stainless steel; in case of contact please rinse immediately.

Warranty

This instrument is designed and built according to the Community Directive CEE 93/42 "*Medical equipments*" and C.B.M. is responsible for safety, reliability and performances provided that possible modifications or repairs are carried out by authorised people, the corresponding electric plant complies with IEC regulations and the equipment is used in compliance with these instruction for use.

C.B.M. guarantees its sterilizer for two (2) years, in case of manufacturing defects or failures of components, not due, according to C.B.M. evaluation, to the operators' wrongful use of the apparatus.

The warranty does not cover parts subject to wear and repairs or replacements carried out by personnel not authorised by C.B.M.

IMPORTANT: The manufacturing company is released by any responsibilities for damages to people or things due to failure to earth the plant.

Details, descriptions and pictures are merely indicative.

C.B.M. S.r.l. reserves the right to carry out ameliorative modifications to its own products without any previous notice.

Disposal Of Electric And Electronic Refuse

Waste from Electrical and Electronic Equipment - Directive 2002/96/EC (WEEE)

The symbol (shown above) that appears on this product means that it should not be disposed of together with normal household refuse but should be treated separately.

Electric and electronic equipment may contain substances that are harmful to the environment and to the health of humans and should therefore be disposed of at special waste disposal sites.

There are harsh legal penalties in force against persons found disposing of such equipment in an irresponsible manner.

When the equipment is no longer of any use, it may taken back to the outlet from which you intend to purchase a replacement.

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