User Manual

Emergency Trolley

(EU)MDR2017/745

XH-75071A-LS

Basic UDI-DI:6974580870009E6

I. Product Overview

This emergency cart is a mobile medical device specifically designed for healthcare institutions. It facilitates the rapid and convenient storage and transport of emergency medications, instruments, and other supplies to meet urgent medical treatment needs, thereby securing valuable treatment time for patients.

II. Structural Components

1. Main Frame: JEMP assembly line product primarily constructed from steel and ABS engineering plastic; aluminum alloy triangular pillars provide four-column load-bearing support.
2. Upper Body: ABS curved base formed via injection molding with side handrails; recessed design prevents item slippage. Work surface features stainless steel guardrails and tempered glass overlay for scratch resistance and easy cleaning.
3. Front Section: The central locking system is foldable. The vehicle body features a total of five drawers, each equipped with adjustable dividers that can be installed or removed based on actual usage requirements.
4. Left side: Defibrillator platform with reinforced straps for securing the device; retractable auxiliary work surface for writing, which slides in for space-saving when not in use; Removable file box (also called document box) on the left serves as a workspace; documents can be stored inside after use;
5. Right side of the cart: Sharps container and dual-color waste bins. The right side is the contamination zone. Discarded needles go into the sharps container for single-use disposal. Two waste bins facilitate waste sorting. Inserting a plastic bag inside makes usage more convenient;
6. Rear of cart: Defibrillation pad (also called compression pad) placed under the patient's

back during resuscitation to ensure stable and safe chest compressions. Oxygen cylinder holder accommodates a 4-liter oxygen tank;

7. Undercarriage: Premium swivel casters with silent rollers, two equipped with brakes (to secure the cart and prevent unintended movement).

III. Usage Procedures

1. Preparation

- Inspect all cart components for integrity. Verify smooth drawer operation, functional castor brakes, and secure IV pole mounting.
- Organize emergency medications, instruments and supplies into designated drawers and platforms based on clinical needs, ensuring items are neatly arranged and easily accessible.

2. Movement and Stabilization

- To move the cart, release the castor brakes and utilize the swivel wheels for effortless maneuvering. Upon reaching the destination, engage the castor brakes to secure the cart.

3. Item Access

- To access items in drawers, pull the corresponding red handle to open the drawer. Return the drawer to its original position after retrieving items.
- The top platform allows direct placement or operation of relevant instruments and equipment.
- Adjust the IV pole height as needed to hang infusion bags for treatment.

4. Waste Disposal

- Dispose of medical waste by sorting it into the corresponding color-coded waste bins according to regulations. Clean and disinfect waste bins promptly after use.

IV. Maintenance and Care

1. Cleaning: Regularly clean the ambulance cart's surfaces, drawers, waste bins and other components using mild detergent and a damp cloth to ensure hygiene and prevent cross-contamination. Dry thoroughly with a clean cloth afterward to avoid rust caused by moisture residue.

2. Inspection: Regularly inspect drawer slides, casters, IV stands, and other components. Apply lubricant if slides are stiff. Repair or replace worn or unresponsive casters promptly. Secure any loose IV stands.

3. Storage: When the ambulance is not in use for extended periods, store it in a dry, well-ventilated area away from direct sunlight to prevent environmental damage.

V. Precautions

1. Regularly inspect medications and instruments inside the ambulance to ensure they remain within their expiration dates and in good working order. Replenish depleted supplies promptly.

2. When moving the ambulance, be mindful of the surrounding environment to avoid collisions with walls, equipment, or personnel.

3. During operation, gently pull and push drawers, auxiliary workstations, and other components to prevent damage from excessive force.

Unauthorized personnel must not disassemble or modify ambulance components to avoid compromising normal operation and performance.

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Zhangjiagang Xiehe Medical Apparatus & Instruments Co., Ltd No.7th, Middle Xinzha Road,
SUNGO Europe B.V.