



User Manual



Stair Stretcher YXH-5C

(EU) MDR 2017/745
Basic UDI-DI: 6974580870001DN



Accessories:

Three safety belts.

Technical data:

- 1.Product size(L×W×H):84×60×98cm, Folding size(L×W×H):98×60×20cm
- 2.Cartoon size(1pc):100×61×20cm
- 3.N.W.(1pc): 14kg G.W.:15.5kg
- 4.Load bearing: ≤159kg

Description:

1. It is mainly made of high-strength aluminum alloy materials .
2. It has four wheels and will be easy for moving on flat ground.
3. It is characterized by its being light-weight, small-sized, easy-carried, use-safely and easy for disinfection and cleaning.

Intended Use:

It is designed for transferring people with limited mobility on stairs and flat ground.

Attention:

This device shall only be operated by fully trained persons,untrained persons are prohibited from operating this device.

Maximum load capacity: 159 kg. Do NOT exceed the weight limit.

The safety belt must always be fastened during use.

Inspect the equipment before every use. Check for loose or damaged components.

Do NOT use the product if any defect is found.

Maintenance:

Keep clean as routine (include disinfection). Compatible disinfectants: 75% ethanol, chlorine-based disinfectant (follow manufacturer's concentration instruction). Do NOT use corrosive chemicals.



Zhangjiagang Xiehe Medical Apparatus & Instruments Co., Ltd
No.7th, Middle Xinzha Road, Zhashang Industrial Zone,Yangshe Town, Zhangjiagang City, Jiangsu
215600, China.
SRN: CN-MF-000008449



Quirumed S.L.U. CIF B97267405 C/ Corretger 117, Pq. Emp.Táctica, Paterna, Valencia (Spain).



SUNGO Europe B.V.
Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands.
SRN: NL-AR-000000247