User Manual

Emergency Trolley (EU)MDR2017/745

XH-85037 Basic UDI-DI:6974580870009E6

I. Product Overview

This emergency cart is a mobile medical device specifically designed for healthcare institutions. It facilitates the rapid and convenient storage and transport of emergency medications, instruments, anesthetic drugs and other supplies to meet urgent medical treatment needs, thereby securing valuable treatment time for patients.

II. Structural Components

It is mainly composed of steel and ABS engineering plastic structure. Aluminum alloy triangular and four-column load-bearing.

2. Upper part of the vehicle body: The platform is in a flying swallow style, with ABS guardrails formed in one piece. The gentle slope design makes it easy to clean and prevents items from slipping. The height of the guardrail is 65mm. The dimensions used inside the platform are 520*430mm. The platform is equipped with transparent soft glass. The left and right handrails of the platform are double-layered, and the concave handprints on the handles are more beautiful and comfortable.

3. Left side of the vehicle body: Embedded telescopic design sub-worktable, side-hanging container box;

4. On the right side of the vehicle body: A rotating circular sharp tool box, a flip-top color-separated trash can, and a two-color trash can for garbage classification.

5. Top of the vehicle body: Adjustable double-row ten-compartment transparent container box (convenient for observation and quick medicine retrieval);

6. Front of the vehicle body: The central control lock is rotatable, equipped with five layers of drawers of different sizes and three-fold soft guide rails. Inside the drawers, there are 3*3 partition plates that can be freely separated. The drawer handle is of the dovetail type, thickened and enlarged, with an ergonomic inclined design.

III. Usage Procedures

1. Preparation

◦ Inspect all cart components for integrity. Verify smooth drawer operation and functional castor brakes.

◦ Organize emergency medications, instruments anesthetic drugs and supplies into designated drawers and platforms based on clinical needs, ensuring items are neatly arranged and easily accessible.

2. Movement and Stabilization

◦ To move the cart, release the castor brakes and utilize the swivel wheels for effortless maneuvering. Upon reaching the destination, engage the castor brakes to secure the cart.

3. Item Access

◦ To access items in drawers, pull the corresponding red handle to open the drawer. Return the drawer to its original position after retrieving items.

◦ The top platform allows direct placement or operation of relevant instruments and equipment.

4. Waste Disposal

◦ Dispose of medical waste by sorting it into the corresponding color-coded waste bins according to regulations. Clean and disinfect waste bins promptly after use.

IV. Maintenance and Care

1. Cleaning: Regularly clean the ambulance cart's surfaces, drawers, waste bins and other components using mild detergent and a damp cloth to ensure hygiene and prevent cross-contamination. Dry thoroughly with a clean cloth afterward to avoid rust caused by moisture residue.

2. Inspection: Regularly inspect drawer slides, casters and other components. Apply lubricant if slides are stiff. Repair or replace worn or unresponsive casters promptly.

3. Storage: When the ambulance is not in use for extended periods, store it in a dry, well-ventilated area away from direct sunlight to prevent environmental damage.

V. Precautions

1. Regularly inspect medications and instruments inside the ambulance to ensure they remain within their expiration dates and in good working order. Replenish depleted supplies promptly.

2. When moving the ambulance, be mindful of the surrounding environment to avoid collisions with walls, equipment, or personnel.

3. During operation, gently pull and push drawers, auxiliary workstations, and other components to prevent damage from excessive force.

Unauthorized personnel must not disassemble or modify ambulance components to avoid compromising normal operation and performance.

How to install two anesthetic storage cabinets:

1. Each trolley is equipped with two cabinets for storing anesthetics, Please be sure to distinguish between the large and small holes: the large holes are circled in red, and the small holes are circled in blue. The tubes should be inserted into these four holes.

2. Pick up the tube and insert the **the threaded section end(as the below picture)** into the four holes marked in the previous step.

3. **Install Cabinet A first, then install Cabinet B.**
4. Once installation is complete, use the included screws to secure the top.

Slide the washer onto the screw, screw the screw into the hole at the top (please use the tool provided by the manufacturer).

5. Insert the assembled unit into the trolley body.

Unscrew the plastic cap from the connector, insert the cap into the tube first, then insert the tube into the connector, lower the cap, and tighten it by hand.

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EC REP SUNGO Europe B.V.