Camilla funeraria plegable

Modelo: YXH-1H (Referencia: 1238-YXH1H)

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Model: YXH-1H (Ref: 1238-YXH1H)

Technical data:

Material: High strength aluminum alloy and PVC.
Product size: 205×53×35cm.
Folding size: 108×53×10cm.
Packing size: 112×54×12cm (1pc/carton).
Load Capacity: 159kg.
Net weight: 9kg. Gross weight: 10kg.

Operation Procedure:

- Remove the stretcher and place it flat on a smooth surface.
- Press down on the button (as shown in the picture ①) to unfold the wheels and legs and position them on the ground.
- Open the top cover of the stretcher and lift the body onto it. Secure with safety belts , ensure the weight of the body is within the stretcher's load capacity.
- Close the top cover of the stretcher. Have someones stand at the two ends or at the four corners of the stretcher to lift it steadily and transport.
- For short-distance, flat-ground transportation, raise the foot pedals(as shown in the picture ②) on the legs of the stretcher. Keep the small wheels at the rear end(as shown in the picture ③) on the ground. One person can smoothly push it forward. Maintain the stability of the stretcher throughout the process to avoid the body from tilting.

Applied for:

The funeral folding stretcher YXH-1H is mainly used in the hospitals and funeral houses to transport dead bodies.

Attention:

- When carrying bodies with this stretcher, be sure fasten the safety belts and not exceed the load capacity (159kg).
- In general operation, be careful not to out torn the stretcher surface.

Maintenance：

- Keep clean as routine (incl. sterilization).
- Often check if any parts loose or not.

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Zhangjiagang Xiehe Medical Apparatus & Instruments Co., Ltd