User Manual

Stair Stretcher **(EU) MDR 2017/745**
YXH-5E **Basic UDI-DI:6974580870001DN**



Accessories:
Two safety belts.

Technical data:
1. Product size(L×W×H):92×52× 105cm, Folding size(L×W×H):114×52×26cm
2. Carton size(1pc):115×54×27cm
3. N.W.(1pc): 12kg G.W.:13.5kg
4. Load bearing: ≤159kg

Description:
1. Single-person operation is permitted on flat ground and when descending stairs.Two operators are mandatory for ascending stairs.
2. It has four wheels and will be easy for moving on flat ground.
3. The stretcher approves the patient's safety during transport process with two pieces of safety belts.
4. It is made of aluminum alloy materials. It is characterized by its being light-weight, small-sized, easy-carried, use -safely and easy for disinfection and cleaning.

Intended Use：

It is designed for transferring people with limited mobility on stairs and flat ground.

Attention:
This device shall only be operated by fully trained persons,untrained persons are prohibited from operating this device.
Maximum load capacity: 159 kg. Do NOT exceed the weight limit.
The safety belt must always be fastened during use.
Inspect the equipment before every use. Check for loose or damaged components.
Do NOT use the product if any defect is found.

Maintenance：
Keep clean as routine (include disinfection). Compatible disinfectants: 75% ethanol, chlorine-based disinfectant (follow manufacturer’s concentration instruction). Do NOT use corrosive chemicals.

Zhangjiagang Xiehe Medical Apparatus & Instruments Co., Ltd
No.7th, Middle Xinzha Road, Zhashang Industrial Zone,Yangshe Town, Zhangjiagang City, Jiangsu 215600, China.
SRN: CN-MF-000008449

Quirumed S.L.U. CIF B97267405 C/ Corretger 117, Pq. Emp. Táctica, Paterna, Valencia (Spain).

EC REP
SUNGO Europe B.V.
Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands.
SRN: NL-AR-000000247

Document Revision: Rev.02
Revision date: 2026-10-09