

Patient protection eyewear "ALLROUND", color white

Product information:

Item no.:	600-ALLROUND-02
Frame color:	weiß
Product weight:	0 g



Certified protection specifications for EN 207

Wavelengths	Protection levels
• 180-315	D LB10 + IR LB5 (OD10+)
• >315-1400	D LB7 + I LB9 + R LB8 + M LB10 (OD10+)
• >1400-11000	DI LB5 (OD5+)
	PF S CE

Properties:

Patient goggles for all laser wavelengths made of stainless steel with opaque, non-reflective surface material for maximum safety. Very high wearing comfort due to variably adjustable headband.

The additional rim rubbers ensure a firm and comfortable hold.

- Application: Patient goggles (eyeball goggles) for all laser wavelengths.
- Material: stainless steel
- Color: white
- with turquoise silicon rim for improved wearing comfort plus 1 spare pair enclosed

Certified according to **EN 207**.

Included in delivery: ALLROUND-02 patient goggles, eye rings (4x), cord, case, instructions for use