

Manual Bed

SK057-1

User's Manual



Picture for reference only

Contents

1. Manual introduction.....	2
2. Features.....	2
3. Technical Parameters.....	2
4. Warning.....	5
5. Maintenance.....	6
5. Storage.....	6
6. Transport.....	6

Note: The appearance of products be described in this manual is a little difference`, this is not effect product performance. Please understand if take you some inconvenience on the operation.

1. Manual Introduction

Please read the instructions carefully as they provide important information on SAFETY DURING USE AND MAINTENANCE PROCEDURES descriptions and images contained in the publication are not binding; we have the right to modify any information for enhancement purposes without any obligation to update this manual. Graphics and images contained in this manual are only examples and may improve products.

2. Features

This product is mainly used for hospitals, nursing homes and public spaces

3. Technical parameters

- **Stainless steel headboard and footboard, blue density board**
- **Mesh bed surface**
- **SPR hardened steel rectangular tube structure**
- **Dimensions: 2050 x 910 x 500mm**
- **Drainage bag hooks**
- **Stainless steel IV poles and stands**
- **Backrest adjustment: 0-75° (±10°)**

4. Warning

- ⚠ DO NOT ALLOW SMALL CHILDREN ON OR NEAR BED DURING OPERATION.
DO NOT ALLOW THIS DEVICE TO BE OPERATED BY SMALL CHILDREN.
- ⚠ WHEN OPERATING THE FUNCTION OF THE BED, ALWAYS ENSURE THAT THE INDIVIDUAL CONFINED TO THE BED IS POSITIONED PROPERLY WITHIN THE CONFINES OF THE BED.
DO NOT LET ANY EXTREMITIES PROTRUDE OVER THE SIDE OR BETWEEN THE BED RAILS WHEN PERFORMING ANY FUNCTIONS.
- ⚠ WARNING/CAUTION LABELS APPLIED TO THE BED OUTLINE HAZARDS OR UNSAFE PRACTICES THAT COULD RESULT IN PERSONAL INJURY AND/OR PROPERTY DAMAGE.
- ⚠ WHEN USING NASAL OR MASKED TYPE ADMINISTERING EQUIPMENT, OXYGEN OR AIR TUBING MUST BE ROUTED AND SECURED PROPERLY TO ENSURE THAT TUBING DOES NOT BECOME ENTANGLED AND EVENTUALLY SEVERED DURING NORMAL OPERATION OF BED.
- ⚠ KEEP ALL MOVING PARTS FREE OF OBSTRUCTIONS (I.E. BLANKETS/SHEETS, HEATING BLANKETS/PADS, TUBING, WIRING, AND OTHER TYPES OF PRODUCTS)

5. Maintenance

CHECKLIST

- 5.1 Periodically clean this bed (including disinfection).
- 5.2 Check frequently for loose or damaged parts.