

Luxury Transfusion Chair

SKE004-1

User's Manual



Picture for reference only

1. Manual Introduction

This manual is a part of the article to which it is referred. Please read the instructions carefully as they provide important information on SAFETY DURING USAGE AND MAINTENANCE PROCEDURES descriptions and images contained in the publication are not binding; we has the right to modify any information for enhancement purposes without any obligation to update this manual. Graphics and images contained in this manual are only examples and may improvement of products.

2. Applicable Scope

It is mainly used for hospital, nursing home, home care.

3. Features

- Size: 710*800*102mm
- Powder coated steel structure.
- PU foam seat and back
- Painted wood arm



4. Warning and Caution

- ⚠ DO NOT ALLOW SMALL CHILDREN ON OR NEAR BED DURING OPERATION.
- ⚠ DO NOT ALLOW THIS DEVICE TO BE OPERATED BY SMALL CHILDREN.
- ⚠ WHEN OPERATING THE HI-LO, KNEE, TRENDLENBERG, REVERSE TRENDLENBERG OR BACK FUNCTION OF THE BED, ALWAYS ENSURE THAT THE INDIVIDUAL CONFINED TO THE BED IS POSITIONED PROPERLY WITHIN THE CONFINES OF THE BED. DO NOT LET ANY EXTREMITIES PROTRUDE OVER THE SIDE OR BETWEEN THE BED RAILS WHEN PERFORMING ANY FUNCTIONS.
- ⚠ DO NOT USE UNAUTHORIZED SIDE RAILS.
- ⚠ WARNING/CAUTION LABELS APPLIED TO THE BED OUTLINE HAZARDS OR UNSAFE PRACTICES THAT COULD RESULT IN PERSONAL INJURY AND/OR PROPERTY DAMAGE.
- ⚠ WHEN USING NASAL OR MASKED TYPE ADMINISTERING EQUIPMENT, OXYGEN OR AIR

TUBING MUST BE ROUTED AND SECURED PROPERLY TO ENSURE THAT TUBING DOES NOT BECOME ENTANGLED AND EVENTUALLY SEVERED DURING NORMAL OPERATION OF BED.

⚠ BEFORE OPERATING, INSURE THAT THE BED FRAME RIVETS ARE PROPERLY SEATED IN THE BED END HOOKS: OTHERWISE INJURY OR DAMAGE MAY RESULT.

⚠ KEEP ALL MOVING PARTS FREE OF OBSTRUCTIONS (I.E. BLANKETS/SHEETS, HEATING BLANKETS/PADS, TUBING, WIRING, AND OTHER TYPES OF PRODUCTS)

5. Maintenance

5.1 Keep clean as routine (Incl. disinfect).

5.2 Check if any parts loose or not as usual.

6. Storage


Store this hospital examination couch in the place of damp proof and non-corrosion environment

7. Transport

7.1 Be transferred by common transporting tools: ship, plane, train, truck.

7.2 Hospital bed can be manually handled.

7.3 No throwing or dropping when transshipment.



.....Because we care