

ABS Bedside Cabinet

SKS002

User's Manual



Picture for reference only


Specifications:

- Size: 480*480*760mm
- Frame: ABS Injection.
- Cabinet top, drawer, tray, layer board, door and the base are made of ABS
- Towel hangers on the 2 sides of cabinet, that can be adjustable
- Dinning board is between the drawer and tabletop
- Thermos , washbasin and other commodity can be put in it

Front side



Adjustable plate inside the cabinet



Staff can be put into cabinet door



Towel rack



Retractable dinning table



Dinning table, drawer, cabinet open





MODEL:SKS002 | REF.:319-SKS002

MESA HOSPITALARIA/HOSPITAL TABLE / TABLE DE CHAMBRE D'HOPITAL/
TAVOLINO PER DEGENZA/HOSPITAL DE MESA /KRANKENHAUS TISCH



Jiangsu Saikang Medical Equipment Co.,Ltd

No. 35 Lehong Road, Modern Agriculture
Demonstration Park, Zhangjiagang City,
Jiangsu Province, China



SUNGO Europe B.V.

Fascinatio Boulevard 522, Unit 1.7,
2909VA Capelle aan den IJssel,
The Netherlands



QuirumedS.L.U. CIF B97267405

C/Corretger 117 Pq.Emp.Tactica,Paterna,
46980, Valencia(Spain)



LOT:HW25033



2025-xx-xx